



# dermalight® 800

The versatile full-body UV-irradiation unit



## Therapy in standing, lying or sitting position

The irradiation unit of the dermalight® 800 is continuously slewable by 180° and enables therapy in standing, lying and sitting position. But also local treatment, f. ex. in the neck area, can be carried out without problems. The continuous height adjustment grants always the correct treatment distance. The moveable tripod of the unit grants an easy movement to the treatment place and a space-saving storage.

## Exact dosage with the timer

The treatment times can be dosed exactly with the separate electronic timer, which optimizes the therapy results. The unit is switched off automatically after having reached the adjusted therapy time. In case of interruptions the remaining time is stored. Thanks to the high UV-intensity, the treatment times are reduced, which reduces the formation of erythema and the cumulative dose.



## Indications

UV-B/311nm

- Psoriasis
- Vitiligo

UV-A

- Neurodermatitis
- Acne

## Technical Data

Dimensions (LxWxH):  
800 x 800 x 1900 mm

Weight of the unit:  
approx. 45 kg

Power supply:  
230V / 50Hz / 20VA

## Lamp equipment

UV-B/311nm narrowband  
Art.-no. 500001

UV-A  
Art.-no. 510001



CE 0123



## Special characteristics

- Therapy possible in sitting, lying or standing position
- Secure and easy height adjustment with the pneumatic spring
- Short therapy times
- Moveable and space-saving tripod
- Exact dosage thanks to the separate LCD-timer
- Low power consumption
- With line filter for the compensation of electromagnetic compatibility
- 2 years warranty

## Dr. K. Höhle Medizintechnik GmbH

Viktor-Frankl-Straße 20  
86916 Kaufering, Germany

Phone 0049-8191-30555-0  
Fax 0049-8191-30555-50

e-mail [medizin@drhoenle.de](mailto:medizin@drhoenle.de)  
[www.drhoenle.de](http://www.drhoenle.de)